



PRO-Craft

DENTAL LABORATORY, INC.

Creating Life Changing Experiences

25791 Jefferson Ave. Murrieta, Ca. 92562

(951) 304-2866 • 1-877-4-The-Lab • 1-877-484-3522

pro-craft.com • Info@pro-craft.com

Dr. Name _____ Account # _____

Address _____ City _____ State _____ Zip Code _____

Phone _____ Prep Date _____ Deliver By _____

Patients Name _____

PLEASE SEND: BOXES RX FORMS SHIPPING LABEL/BAGS

Teeth

Premium Teeth (per card) - Additional Fee Anterior Posterior

Basic Teeth (per tooth) Anterior Posterior

Shade _____ Mould _____ Tissue Shade _____

Papillameter _____ Alameter _____

Processing

- Custom Tray
- Model Mount & Registration
- Bite Rim Full Denture
- Reset Denture
- Bite Rim Partial
- 1 Saddle
- 2 Saddles
- 3 Saddles

Complete

- Full Denture Immediate
- Vitallium Unilateral
- Chrome Partial
- Valplast Partial
- Chrome Unilateral
- Vitallium Frame
- Stay Plate
- Valplast Unilateral
- Valplast Metal Frame/Combo
- Acrylic Unilateral

Full Denture Try In w/Teeth

- Full Denture Immediate Try in w/Teeth
- Full Denture Try in w/Teeth
- Full Premium Denture Try in w/Teeth
- Screw Retained Acrylic Hybrid Try In w/Teeth

Full Denture Finish

- Finish Full Denture Immediate
- Screw Retained Acrylic Hybrid Finish
- Finish Full Denture
- Finish Full Premium Denture

Partial Denture Try In w/Teeth

- Chrome Partial Try In w/Teeth
- Unilateral Chrome Try In w/Teeth
- Vitallium Partial Try In w/Teeth
- Unilateral Vitallium Try In w/Teeth
- Valplast Partial Try In w/Teeth
- Unilateral Valplast Try In w/Teeth
- Valplast Metal Frame/Combo Try In w/Teeth
- Unilateral Acrylic Try In w/Teeth
- Stay Plate Try In w/Teeth

Partial Denture Finish

- Chrome Partial Finish
- Unilateral Chrome Finish
- Unilateral Valplast Finish
- Vitallium Partial Finish
- Unilateral Vitallium Finish
- Unilateral Acrylic Finish
- Valplast Metal Frame/Combo Finish
- Valplast Partial Finish
- Stay Plate Finish

Custom Extras

- Ball Clasp
- Name In Appliance
- Hader Bar
- Clear Clasp
- PRO Sil Gasket
- ERA Attachment
- Y Clasp
- Verification Jig
- Bredent Attachment
- Wire Clasp
- Surgical Guide
- Locator Attachment
- Wire Mesh
- Surgical Guide w/Drill Sleeve
- Wire Braided
- Wire Reinforcement
- Wire Mesh Arch
- Duplicate Models

Reline/Repair Denture or Partial

- Soft Reline
- Reline
- Repair
- Wire Mesh
- Rebase
- Weld
- Valplast Repair

Snoring Appliance

- EMA SleepBetter

Lab Time: 3 Days

Impact Zone™ Nightguards


- Relaxer - Anti Clenching
- Dual Laminate
- Soft
- Talon Plus
- Bleach Tray
- Talon Thermoplastic
- Hard Acrylic
- Clear Retainer

Lab Time: 3 Days

Monster™ Sports Mouth Guard

- Junior
- Pro
- Elite
- Pro-Elite

Lab Time: 3 Days



LAB USE ONLY

- Disinfected by: _____
- Impression Plastic Metal Double Bite Study Model
- Bite Yes No
- Poured Models Working Model Oposing
- Number of Old Crowns _____
- Articulator _____
- Sep Dies _____
- Shade Tab _____
- Soft Tissue _____
- Photos _____
- Locator Attachment _____
- Partial _____
- Denture _____
- Digital Model _____
- Alloy _____
- Abutment _____
- Analog _____
- Impression Coping _____
- Screw _____
- Locator Housing _____
- Locator Analog _____
- Custom Ti-Base _____
- Analog/Abut _____
- Repositional Analog _____

- QC Digital File
- Please Call
- Metal Try-In
- Bisque Bake Try-In
- Dr. Die Trim

Note: Please send a patient approved study model on all work involving anterior teeth.

- Remake
- Adjustmet

Shade _____

Stump Shade _____

Implant Information _____

Brand _____

Size _____

Signature _____ D.D.S. License No. _____

Customer agrees to company policy as stated on reverse

Porcelain Fused to Metal Crowns

- Creation Porcelain Fused to Metal
- Non-Precious Alloy
- Semi-Precious Alloy
- Premium Ceramic Alloy 86%

Lab Time: 7 Days

All Ceramic

- All Z[®] Solid Zirconia
- All Z[®] Plus Translucent Solid Zirconia
- All Z[®] Plus Inlay / Onlay
- All Z[®] Zirconia Hybrid
- e.max Crown
- e.max Veneers
- e.max Inlay / Onlay
- Empress Crown
- Empress Esthetic Veneer
- Pro Z Layered Zirconia
- Enamic Vita Composite Crown
- Enamic Vita Composite Inlay/Onlay

Lab Time: 7 Days

Diagnostics

- Triad Diagnostic Wax-up
- Pro-Temps
- Essex Removable Temps
- Snap on Smile
- Digital Diagnostic

Lab Time: 5 Days

Full Metal C&B

- Gold Crown / Inlay / Onlay
- 52% Alloy
- 75% Alloy
- 2% Alloy

Lab Time: 5 Days

Implants

- Screw Retained Crown
- Titanium Custom Abutment
- Zirconia Custom Abutment
- Prep Abutment
- Abutment Placement Stent
- Implant Parts-Manufacture

Lab Time: 7 Days Abutment & Crown

Specialties

- Porcelain 180° Margin
- Rest
- PD Attach
- Soft Tissue
- Seperate Cast Post
- ERA
- Bredent Attach
- Custom Shade
- Adapt to Crown
- Pink Tissue
- Stress Breaking
- Wing
- Hader Ball

If No Occlusal Clearance

- Metal Occlusion
- Spot Opposing
- Reduction Coping
- Non-Vital Prep

Occlusal Staining

- Light
- Medium
- Dark
- None

Pontic Design



Ridge Relief Yes No _____ mm.